

OXYMED




Oxymed P2 Hybrid

Portable Oxygen Concentrator

Flow: Pulse & Continuous


Simplicity Meets Comfort

The 2.8-inch vivid LCD display and simple controls make operation effortless, ensuring comfort and confidence every step of the way.



Seamless Mobility with Roller Cart

Designed for smooth mobility and convenience, the roller cart ensures you stay active and comfortable, wherever life takes you.



Highest Flow Ever

Offers pulse dose delivery with up to 10 settings, accommodating patients with needs of up to 6 liters per minute, and continuous flow delivery of up to 3 liters per minute.

Pulse Flow Setting	Equivalent Flow (LPM)	Battery Backup
Setting 1	0.6	9.50 hours
Setting 2	1.2	8.96 hours
Setting 3	1.8	8.62 hours
Setting 4	2.4	7.04 hours
Setting 5	3.0	6.51 hours
Setting 6	3.6	6.16 hours
Setting 7	4.2	4.40 hours
Setting 8	4.8	3.52 hours
Setting 9	5.4	2.99 hours
Setting 10	6.0	2.50 hours

Continuous Flow Setting	Flow Delivered	Battery Backup
Setting 1	1.0 LPM	6.7 hours
Setting 2	1.5 LPM	4.3 hours
Setting 3	2.0 LPM	3.2 hours
Setting 4	2.5 LPM	2.4 hours
Setting 5	3.0 LPM	2.0 hours

Hybrid Flow Modes

The P2 Hybrid features two oxygen flow modes

Pulse Mode: Recommended for daily activities and mobile use, delivering oxygen only when you inhale, conserving oxygen and extending battery life.

Continuous Mode : Recommended for high oxygen needs, neuromuscular disorders, or CPAP users requiring a steady flow of oxygen.

Disclaimer: Equivalent flow values are indicative and calculated by assuming I:E of 1:2, please consult your doctor before buying the machine




Specifications

Weight	8.0 ± 0.1kg
Size	28.2W x 17.1D x 38.2H cm (includes battery)
Noise Level	49 dB
Flow settings	
Oxygen Flow Output	Pulse: Equivalent to 6 LPM (10 settings) Continuous: 3 LPM (5 settings)
Oxygen Flow Output	Pulse Dose Setting 1 - 200 mL/min (~ 0.6 LPM) Setting 2 - 400 mL/min (~ 1.2 LPM) Setting 3 - 600 mL/min (~ 1.8 LPM) Setting 4 - 800 mL/min (~ 2.4 LPM) Setting 5 - 1000 mL/min (~ 3.0 LPM) Setting 6 - 1200 mL/min (~ 3.6 LPM) Setting 7 - 1400 mL/min (~ 4.2 LPM) Setting 8 - 1600 mL/min (~ 4.8 LPM) Setting 9 - 1800 mL/min (~ 5.4 LPM) Setting 10 - 2000 mL/min (~ 6.0 LPM) ~ signifies the equivalent continuous flow value with I:E of 1:2
	Continuous Flow Setting 1 - 1.0 LPM Setting 2 - 1.5 LPM Setting 3 - 2.0 LPM Setting 4 - 2.5 LPM Setting 5 - 3.0 LPM
Oxygen Concentration	87%-96% at all settings
Warm up Time	5 minutes
Display Screen	2.8 inch

Oxymed P2 Hybrid Portable Oxygen Concentrator

Battery & Power	
Battery Backup	Battery (Included) Upto 9.5 hours
Recharge Time	8 hours
Battery Specs	Battery: Voltage: 14.4VDC Rated capacity: 20Ah
Power Supply	AC Input: 100 to 240VAC 50 to 60 Hz
Warranty & Quality Standards	
Quality Standards	CE Certified
Warranty	3 years (Main unit) 1 year (Battery and sieve beds)
What's in the Box	Oxymed P2 Hybrid AC Adapter Car Charger Battery Nasal Cannula Roller cart Filter set Carry Bag User Manual



Registered Office:
A brand by Medequip Healthcare Solutions Pvt. Ltd
Plot No: A-12, A-13, Apparel Park
Phase-1, Doddaballapur Industrial Area
Bengaluru - 561203, Karnataka

Branch Offices:
Delhi, Telangana, Odisha, Gujarat, Maharashtra,
Uttar Pradesh, Rajasthan, Madhya Pradesh Punjab,
Tamil Nadu, West Bengal, Kerala, Jharkhand,
Chhattisgarh and Bihar

