

BPL Respiratory solutions

Every Breath Matters



Follow Us on



Happier Living Everyday

BPL *Harmony 1*

Auto C-PAP solutions

BPL's Auto CPAP - is a new generation sleep and respiratory therapy platform that offers a systematic solution for Obstructive Sleep Apnoea apart from offering a world class algorithm, which effectively manages the apnoea's and hypopneas. The ergonomic design, with an intuitive interface and a complete range of accessories brings about a brilliant & comfortable experience for users.



Effective Treatment in Obstructive Sleep Apnea with BPL Harmony 1 Auto CPAP solution

Effective Therapy

- Advanced Algorithm
- Data management
7 days of detailed data and 30 days of therapy data
- Complete Range of Comfortable Face Masks
- Quiet Therapy



Easy to use

- Compact and Smart
- Auto On/Off
- Intuitive user interface with an easy to read angled screen

Greater Comfort

- Smaller footprint
- Ramp & Auto Ramp
- Fully Integrated Humidifier
- Intelligent Pressure Release (iPR)



BPL Harmony 1

Quiet therapy

The new advanced flow generator enables the air to exit/flow axially. This leads to less turbulence and a quieter therapy, ensuring that the patients & their bed partners are not disturbed during their sleep.

Enhanced Clinical algorithm

BPL Harmony1's intelligent algorithm maintains the ideal pressure for patients & automatically responds to apnea, snoring and flow limitations thereby keeping the airway open and ensuring comfortable apnea free sleep for the user throughout the night.

Intelligent Pressure Release (iPR)

BPL Harmony 1's intelligent pressure release (iPR) feature, drops the pressure upon expiration by a defined value thereby ensuring a comfortable breathing experience for the user.

Ramp / Auto Ramp

Ramp enables therapy to begin at a minimum pressure thus making the beginning of treatment more comfortable for the patient. The pressure continues to increase with the gradual ramping up to the minimum therapy pressure within 5 - 60 min mins, coinciding with the patients sleep onset.



Auto On/ Off

Convenient for patients as the treatment automatically starts/stops when and if the mask is taken off or put on after a break during therapy, without any external intervention.

Data Recording

BPL Harmony 1 records up to 365 days of therapy data (Pressure, AHI, Leak, AI, HI, etc) and 7 days of breath by breath detailed data, helping patients and their doctor monitor their therapy for better prognosis

Integrated Humidifier

Effective humidification leads to a better and comfortable sleep without the discomfort of a dry or sore throat. An integrated humidifier allows for this benefit without exponentially increasing the Auto CPAP's footprint.

Colour Display

The colour LCD screen on BPL Harmony 1 with a push dial control knob and an easy-to-read menu lets the user use the device seamlessly, providing an unparalleled user experience.



BPL *LifePAP*

BPL LifePap BiPAP Solutions

BPL Life PAP is a non-invasive ventilator (BiPAP) designed for non-dependent patients with obstructive or restrictive respiratory conditions.

Featuring BPL's volume-assurance mode – eVAPS (Efficient Volume-Assured Pressure Support) – BPL's LifePap maintains target tidal volume to suit each patient's changing needs. With S, ST and ST with AutoEPAP modes, BPL Life PAP is easy to set up and use, supported by BPL's patient management software.





Effective treatment of OSA with the BPL LIFE PAP (BiPAP) solution

Efficient Volume-Assured Pressure Support (eVAPS)

eVAPS or Efficient Volume Assured Pressure Support when activated, will ensure a specified constant average tidal volume "Vt" is effectively delivered. Life PAP will adjust the Pressure Support (PS) to meet the said Vt targets. During the therapy, IPAP is changed to meet the Pressure Support, adjustment demand.

AutoEPAP

A feature only available with a few, Auto EPAP automatically addresses any upper airway obstruction that the patient may be experiencing during the course of the therapy, ensuring minimal physical interventions by family members or paramedics to address challenges of high PCO₂ due to such an obstruction.

Integrated Humidifier

Effective humidification leads to a better and comfortable therapy without the discomfort of a dry or sore throat. An integrated humidifier allows for this benefit without exponentially increasing the BiPAP's footprint.

Built-In Simplicity

The user interface on BPL LifePap NIV has been designed keeping the patients in mind. Intuitive and simple to use, easy to view with personalized patient comfort settings. The Touch screen allows for a seamless navigational experience, a feature that is not available today in most of the other BiPAP's in the market.

A complete solution

BPL Life PAP ST with Auto EPAP modes offers all essentials operating modes for patients like CPAP, S, ST and ST with Auto EPAP. The in-built Alarms are a very useful feature for hospitals who want a state of the art NIV solutions for their patients.

BPL *LifePAP*

Quiet therapy

The new advanced flow generator enables the air to exit/ flow axially. This leads to less turbulence and a quieter therapy, ensuring that the patients & their bed partners are not disturbed during their sleep.

Advanced algorithm

BPL LifePAP's Advanced Ventilation algorithm works in synchrony with the patient's own breathing pattern, thereby improving the patient's compliance with the therapy.

Auto E-pap

Auto E-pap automatically addresses any upper airway obstruction that the patient may be experiencing during the course of the therapy.

Intelligent Pressure release (iPR)

Intelligent pressure release is designed to make the therapy more comfortable by reducing the treatment pressure in the cycle process and expiration. By applying iPR, the average treatment pressure is lowered than its counterpart in normal operation.



Auto On/ Off

Convenient for patients as the treatment automatically starts/stops when and if the mask is taken off or put on during therapy.

Integrated Humidifier

Effective humidification leads to a better and comfortable sleep without the discomfort of a dry or sore throat.

Simple and intuitive user interface

Simple and intuitive user interface with a personalized patient comfort settings and different modes, ensures a seamless and a worry-free therapy experience for the patient. A Touch Screen further allows for a hassle free navigational experience.

e-VAPS (Volume Assured Pressure Support)

e-VAPS intelligently ensures the required tidal volume for patients by automatically adjusting the required pressure support





LIFEPAP Models & Modes

Models	: 25 ST with Auto EPAP	25ST	25S
eVAPS	: ✓	✓	
Auto EPAP	: ✓		
ST Mode	: ✓	✓	
S Mode	: ✓	✓	✓
CPAP	: ✓	✓	✓

Product Specifications

Parameters	BPL Harmony 1 Auto CPAP	BPL LIFE PAP BiPAP
Display	3.5" Colour LCD Display	5" Colour touch screen display
Flow	>120 L/min	>120L/min
Safety Precautions	Class II Type BF	Class II Type BF
Ingress Protection	IP22	IP22
Air Filter Efficiency	>75% (7-micron dust)	
Pressure		
Setting range	4-20 cmH2O	4 -25 cmH2O
Max Single Fault Steady Pressure	40 cmH2O	40cm H2O
Dynamic Pressure variation	< 1 cmH2O	<1 cmH2O
Pressure Monitoring Accuracy	± (2%*full scale reading + 4%*actual reading)	± (2%*full scale reading + 4%*actual reading)
Humidifier		
Water Capacity	290 ml (MAX Water Level)	290ml (MAX Water Level)
Humidity	>10 mg/L BTPS (Within the set pressure range)	>10mg/L BTPS (Within the set pressure range)
Noise Value		
A weighted sound pressure level	≤28dBA (uncertainty of 2 dBA)	≤28dBA (uncertainty of 2 dBA)
A weighted sound power level	≤36 dBA (uncertainty of 2 dBA)	≤36 dBA (uncertainty of 2 dBA)
Dimension		
Dimension	222*125*122mm (L × W × H)	235 × 148× 125mm (L × W × H)
Weight	1.35 kg	1.55Kg
Power		
AC Power Input	100~240 Vac 50/60 Hz 2.0A Max	100~240Vac, 50/60Hz, 2.0A max
DC Power Input	24 Vdc, 2.5A	24 Vdc, 2.5A
Operating conditions		
Temperature	5°C- 35°C	5°C- 35°C
Relative Humidity	10% - 95% (Non condensing)	10% - 95% (Non condensing)
Atmospheric Pressure	70kPa – 106kPa	70kPa – 106kPa
Storage Conditions		
Temperature	-25°C - 70°C	-25°C - 70°C
Relative Humidity	5% - 95% (Non condensing)	5% - 95% (Non condensing)
Atmospheric Pressure	70kPa – 106kPa	70kPa – 106kPa

BPL Range of Accessories (CPAP & BiPAP)



Standard Air Tube 6 Feet



Standard Air Tube 6 Feet



Heated tube



AC Adaptor



Nasal Mask : Small



Nasal Mask : Medium



Nasal Mask : Large



Air Filter



Full Face Mask : Small



Full Face Mask : Medium



Full Face Mask : Large



CERTIFIED ISO 13485 : 2016 COMPANY

BPL Medical Technologies Private Limited
Regd. Office: 11th KM, Bannerghatta Road,
 Arakere, Bangalore - 560076, India.
 Toll Free: 1800-4252355
 Website: www.bplmedicaltechnologies.com
 For Enquiries: sales.medical@bpl.in
 CIN: U33110KA2012PTC067282



Happier Living Everyday