



# Quietest BiPAP

Sleep Better. Breathe Better

A BiPAP device that is quiet and has an improved monitoring system for a restful night's sleep.



# RESPIRARE VT 60**ST**

Respirare BiPAP is a non-invasive ventilator designed for spontaneously breathing non-dependent patients with respiratory insufficiency.

# BIPAP VT 60ST







If leaks occur, alarms are triggered automatically, accompanied by a beeping sound indicating the presence of a leak.



Gender Selection helps to give a Soft/Standard response to the selected gender accordingly.

# **Technical Specification**

Modes CPAP, S, T, ST, VAPS Operating/ Monitoring Range CPAP: 4-20 cmH2O IPAP: 4-30 cmH2O | EPAP: 4-28 cmH2O VT: 150-1800 ml Ti/ Ti min/ Ti max 0.2-4 sec/ 0.2-3 sec/ 0.4-4 sec

Sensitivity Off, 5-45 mins/100-900 ms Ramp Time/ Rise Time Peak Flow 240 LPM **Respiratory Rate** 4-45 BPM **Gender Selection** Male/ Female **Humidification Level** Off, 1-5 level Display 2.8 in LCD Alarms Leak Alert with Beep 1/7/31 days, 6 months, 1 year of detailed data. Data Management Dimensions (LxWxH)/ Weight 280x140x95 mm/ 1700 g (Device with Humidifier) **Power Adapter** Input Range: 100-240V, 50/60 Hz, 3A 90 Watt

Humidity/ Altitude non-condensing/sea level to 2,591 m, air pressure range 735 hpa to 1013 hpa -20 c to +60 c / 5-95% relative humidity, Storage and Transport Temperature/Humidity non-condensing Maximum Flow: 4L/min (CPAP), 15 L in S, T, ST, VAPS Supplement Oxygen Standard Air Filter Pollen Filter (Washable) Tubing/ Air Outlet Standard (19 mm)/ 22 mm Air Outlet (STD) Sound Pressure Level <28 dBA with uncertainty of 2 dBA

(typical power consumption)

+5 c to +35 c/ 10-95% relative humidity,

Warranty 2 Years **Conforms to Standards** IEC 60601-1-1

IEC 60601-1-2 IEC 80601-2-12 CDSCO

Class 2 (double insulated, Type BF, Ingress Protection IP22)

#### \*Term & Conditions Apply

Operating Temperature/

## **Key Features**

#### i-Mode

This feature will turn the device **On** for you have the mask on without the need to press any buttons and will also turn the device Off once you take your mask Off.

## **Integrated Humidification**

comfortable and can reduce side effects like a dry

### Mobility, Portability

Easy to carry the device from one location to another.

### **Product Code**

Respirare BiPAP with Humidifier VT 60 **ST** 

#### Manufactured by:









